

TRAUMA SCISSORS

SPECIFICATION

■ # CPBJ0082

■ # CPBJ0020

SIZE: 7.5in(190mm)

WEIGHT: 42g

MATERIAL: Stainless

DESCRIPTION

RHINO RESCUE Trauma Scissors—a serrated side blade design that allows for easy cutting of fabrics, bandages, gauze, and more. The matte black surface layer is alkaline oxidized based on the original color of stainless steel, isolated from air, and resistant to acid and rust. Buckle design, easy to carry, suitable for various outdoor, military, home and other scenarios.



PP Handle

PP handle, fine work-manship, feel comfortable.

Sharp and Worn Teeth

Serrated side blade design, strong bite force, not easy to slip off, easy to break.

Round Hole Design

Round hole design, can be used to peel wire or cut fine iron wire.

Buckle Handle

Buckle design, easy to carry.

Fillet Design

Head fillet design to avoid secondary injury to the wound.



EU REP SUNGO Cert GmbH,
Harffstr.47,40591 Düsseldorf, Germany
E-mail: de.rep@sungoglobal.com



EMS RUN GROUP
8-7 Huide Road, Huishan Area, Wuxi, 214000
Jiangsu, CHINA
E-mail: info@rhinorescuestore.com



WUXI EMSRUN TECHNOLOGY CO., LTD.
Room 501, No.8-7 Huide Road, Wuxi City, Jiangsu
Province, China
E-mail: info@emsrun.com

www.rhinorescuestore.com | info@rhinorescuestore.com