

BLEEDING CONTROL KIT




SPECIFICATION

PRODUCT SIZE: 260×260mm/10.2×10.2in PRODUCT WEIGHT: 440g

DESCRIPTION

The trauma kit contains necessary combat-proven first-aid devices for the staunching of blood flow from traumatic hemorrhage wounds in pre-hospital emergency situations. The stop the bleed kit provides rapid access to bleeding control products and saves lives of severe hemorrhage in incidents.

PRODUCT DESCRIPTION



- | | |
|---------------------------|--------|
| Hemostatic Gauze 4" × 1yd | ×1 |
| Emergency Bandage 6" | ×1 |
| Metal Tourniquet | ×1 |
| Compressed Gauze | ×2 |
| Emergency Blanket | ×1 |
| Bandage Scissors | ×1 |
| Nitrile Gloves | ×1pair |
| Instruction Card | ×1 |

POINTS FOR ATTENTION

- Place sterile compress or clean substitute (a sock, T-shirt, towel, etc.) directly on the wound and apply steady, firm hand pressure.
- If compress soaks through, DO NOT REMOVE. Apply additional compresses and continue hand pressure and Elevate the limb above the level of the heart.<98% EFFECTIVE>

